

Gel/PCR DNA Fragments Extraction Maxi Kit

For research use only

Sample	: up to 6 g of agarose gel up to 2 ml of PCR products
Recovery	: up to 95%
Format	: spin column
Operation time	: 40 minutes
Elution volume	: 0.5-1 ml

Geneaid



www.geneaid.com

Introduction

The Gel/PCR DNA Fragments Extraction Maxi Kit was designed to recover or concentrate DNA fragments (200 bp→10 Kb) from agarose gel, PCR, or other enzymatic reactions. Chaotropic salt is used to dissolve agarose gel and denature enzymes. DNA fragments in the chaotropic salt are bound by the glass fiber matrix of the spin column (1). Contaminants are removed with a Wash Buffer (containing ethanol) and the purified DNA fragments are eluted by a low salt Elution Buffer or TE. Salts, enzymes and unincorporated nucleotides can be effectively removed from the reaction mixture without phenol extraction or alcohol precipitation. Typically, recoveries are up to 90% for Gel Extraction and up to 95% for PCR Clean up. The entire procedure can be completed in 40 minutes and the eluted DNA is ready for use in PCR, Fluorescent or Radioactive Sequencing, Restriction Enzyme Digestion, DNA Labeling and Ligation. For users who require a higher recovery from small base pair DNA fragments (50-200 bp) or large base pair DNA fragments (> 8 Kb), see the order information below.

Quality Control

The quality of the Gel/PCR DNA Fragments Extraction Maxi Kit is tested on a lot-to-lot basis by isolating DNA fragments of various sizes from either aqueous solutions or agarose gel. The purified DNA is checked by electrophoresis.

Kit Contents

Name	DM002	DM010	DM025
DF Buffer	25 ml	125 ml	240+50 ml
W1 Buffer	20 ml	90 ml	210 ml
Wash Buffer* (Add Ethanol)	5 ml (20 ml)	25 ml (100 ml)	50 ml + 25 ml (200 ml) (100 ml)
Elution Buffer (10 mM Tris-HCl, pH 8.5 at 25°C)	2 ml	12 ml	30 ml
DM Column	2 pcs	10 pcs	25 pcs

Order Information

Product Name	Package Size	Cat. No.
Gel/PCR DNA Fragments Extraction Kit	100/300 preps	DF100/300
Small DNA Fragments Extraction Kit (optimized for 50-200 bp)	100/300 preps	DF101/301
Large DNA Fragments Extraction Kit (optimized for > 8 Kb)	100/300 preps	DF102/302
96-Well Gel/PCR DNA Extraction Kit	2/4/10 X 96 Wells	DFP02/04/10
Vacuum Manifold (Accessories)	1 SET	ZVF01

*Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for details).

Caution

DF Buffer contains guanidine thiocyanate which is a harmful irritant. During operation, always wear a lab coat, disposable gloves, and protective goggles.

References

(1) Vogelstein, B., and Gillespie, D. (1979) Proc. Natl. Acad. Sci. USA 76, 615.

Additional requirements

centrifuge tubes, absolute ethanol

Gel Extraction Protocol

Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).

Step 1 Gel Dissociation	<ul style="list-style-type: none"> ● Excise the agarose gel slice containing relevant DNA fragments and remove any extra agarose to minimize the size of the gel slice (TAE buffer is recommended for gel formation). ● Transfer up to 6 g of the gel slice to a 15 ml centrifuge tube. ● Add 10 ml of DF Buffer to the sample and mix by vortex. ● Incubate at 55-60°C for 10-15 minutes or until the gel slice has been completely dissolved. During incubation, invert the tube every 2–3 minutes. ● Cool the dissolved sample mixture to room temperature.
Step 2 DNA Binding	<ul style="list-style-type: none"> ● Place the DM Column in a 50 ml centrifuge tube. ● Transfer 12 ml of the sample mixture from the previous step to the DM Column. ● Centrifuge at 14-16,000 x g for 2 minutes. ● Discard the flow-through and place the DM Column back in the 50 ml centrifuge tube (If the sample mixture is more than 12 ml, repeat the DNA Binding Step).
Step 3 Wash	<ul style="list-style-type: none"> ● Add 8 ml of W1 Buffer into the DM Column. ● Centrifuge at 14-16,000 x g for 2 minutes and then discard the flow-through. ● Place the DM Column back in the 50 ml centrifuge tube. ● Add 12 ml of Wash Buffer (ethanol added) into the DM Column and let stand for 3 minutes. ● Centrifuge at 14-16,000 x g for 2 minutes and then discard the flow-through. ● Place the DM Column back in the 50 ml centrifuge tube. ● Centrifuge again at 14-16,000 x g for 5 minutes to dry the column matrix.
Step 4 DNA Elution	<ul style="list-style-type: none"> ● Transfer the dried DM Column to a new 50 ml centrifuge tube. ● Add 0.5-1 ml of Elution Buffer or TE into the center of the column matrix. ● Let stand for 3 minutes or until the Elution Buffer or TE is absorbed by the matrix. ● Centrifuge for 5 minutes at 14-16,000 x g to elute the purified DNA.

PCR Clean Up Protocol

Add absolute ethanol to the Wash Buffer prior to initial use (see the bottle label for volume).

Step 1 Sample Prep.	<ul style="list-style-type: none"> ● Transfer up to 2 ml of a reaction product to a 15 ml centrifuge tube. ● Add 5 volumes of DF Buffer to 1 volume of the sample and mix by vortex.
Step 2 DNA Binding	<ul style="list-style-type: none"> ● Place a DM Column in a 50 ml centrifuge tube. ● Transfer the sample mixture from step 1 to the DM Column and centrifuge at 14-16,000 x g for 2 minutes. ● Discard the flow-through and place the DM Column back in the 50 ml centrifuge tube.
Step 3 Wash	<ul style="list-style-type: none"> ● Add 12 ml of Wash Buffer (ethanol added) into the center of the DM Column and let stand for 3 minutes. ● Centrifuge at 14-16,000 x g for 2 minutes. ● Discard the flow-through and place the DM Column back in the 50 ml centrifuge tube. ● Centrifuge again for 5 minutes at 14-16,000 x g to dry the column matrix.
Step 4 DNA Elution	<ul style="list-style-type: none"> ● Transfer the dried DM Column to a new 50 ml centrifuge tube. ● Add 0.5-1 ml of Elution Buffer or TE into the center of the column matrix. ● Let stand for 3 minutes or until the Elution Buffer or TE is completely absorbed by the matrix. ● Centrifuge for 5 minutes at 14-16,000 x g to elute the purified DNA.

Troubleshooting

Problem	Possible Reasons/Solution
Low Yield	<p>Gel slice did not dissolve completely</p> <ul style="list-style-type: none"> ● The Gel slice was too big. If using more than 6 g of gel slice, separate it into multiple tubes. ● Raise the incubation temperature to 60°C and extend the incubation time. <p>Incorrect DNA Elution Step</p> <ul style="list-style-type: none"> ● Ensure that the Elution Buffer is completely absorbed after being added to the center of the DM Column. <p>Incomplete DNA Elution</p> <ul style="list-style-type: none"> ● If the DNA fragments are larger than 10 Kb, use preheated Elution Buffer (60-70°C) in the Elution Step to improve the elution efficiency.
Eluted DNA doesn't perform well in downstream applications.	<p>Residual ethanol contamination</p> <ul style="list-style-type: none"> ● Following the Wash Step, dry the DM Column with additional centrifugation at 14-16,000 x g for 5 minutes or incubate at 60°C for 5 minutes. <p>DNA was denatured (a smaller band appeared on gel analysis)</p> <ul style="list-style-type: none"> ● Incubate the eluted DNA at 95°C for 2 minutes, and then cool down slowly to re-anneal the denatured DNA.
Low A260/A230	<ul style="list-style-type: none"> ● In the wash step, repeat the 12 ml of Wash Buffer addition and let stand for 2 minutes.