



Explore our innovative technologies for *YOUR* research

RAT ALBUMIN ENZYME IMMUNOASSAY KIT

catalogue # A05102


96 wells

TABLE OF CONTENTS

Presentation	2
Precautions for use	2
Principle of the assay	2
Materials and equipment required	3
Sample collection & preparation	3
Reagent preparation	4
Assay procedure	4
Plate preparation	4
Distribution of reagents and samples	4
Pipetting the reagents	4
Incubating the plate	5
Developing and reading the plate	5
Data analysis	6
Typical data	6
Example data	6
Acceptable range	7
Assay validation and characteristics	8
Assay trouble shooting	8
Bibliography	8



Explore our innovative technologies for **YOUR** research



U.S. patent # 50 47 330
European patent # 89 139 552

**THE RAT ALBUMIN
ENZYME IMMUNOASSAY HAS
BEEN DEVELOPED AND
VALIDATED BY SPI-BIO.**

*For research laboratory use only.
Not for human diagnostic use.*



Société de Pharmacologie et d'Immunologie - BIO

Parc d'Activités du Pas du Lac – Bertin Group
10 bis avenue Ampère
F-78180 – Montigny Le Bretonneux

FRANCE

☎: 33 (0)1 39 30 62 60

📠: 33 (0)1 39 30 62 99

E-Mail: sales@spibio.com

Web: www.spibio.com

February 06



Explore our innovative technologies for **YOUR** research

RAT ALBUMIN EIA KIT

96 wells

Storage : -20°C

Expiry date: stated on the package

This kit contains:

- ☞ A covered 96 well plate, pre-coated with mouse anti-rabbit IgG, ready to use after thawing
- ☞ One vial of Rat albumin tracer, lyophilised
- ☞ Two vials of Rat albumin standard, lyophilised
- ☞ One vial of Rat albumin antiserum, lyophilised
- ☞ Two vials of Quality Control sample, lyophilised
- ☞ One vial of EIA buffer, lyophilised
- ☞ One vial of concentrated Wash buffer, liquid
- ☞ One vial of tween 20, liquid
- ☞ Two vials of Ellman's reagent, lyophilised
- ☞ One instruction booklet
- ☞ One template sheet
- ☞ One well cover sheet

Each kit contains sufficient reagents for 96 wells. This allows for the construction of one standard curve in duplicate and the assay of 33 samples in duplicate.

PRECAUTIONS FOR USE

Users are recommended to read all instructions for use before starting work.

Each time a new pipet tip is used, aspirate a sample or reagent and dispense it back into the same vessel. Repeat this operation two or three times before distribution.

For research laboratory use only.

Not for human diagnostic use.

Do not pipet liquids by mouth.

Do not use kit components beyond the expiration date.

Do not eat, drink or smoke in area in which kit reagents are handled.

Avoid splashing.

The total amount of reagents contains less than 100 µg of sodium azide. Flush the drains thoroughly to prevent the production of explosive metal azides.

PRINCIPLE OF THE ASSAY

This Enzyme Immunoassay (EIA) is based on the competition between unlabelled rat albumin and acetylcholinesterase (AChE) linked to rat albumin (tracer) for limited specific rabbit anti-rat albumin antiserum sites.


The complex rabbit antiserum-rat albumin (free albumin or tracer) binds to the mouse monoclonal anti-rabbit antibody that is attached to the well.

The plate is then washed and Ellman's Reagent (enzymatic substrate for AChE and chromogen) is added to the wells.

The AChE tracer acts on the Ellman's Reagent to form a yellow compound.

The intensity of the colour, which is determined by spectrophotometry, is proportional to the amount of tracer bound to the well and is inversely proportional to the amount of free rat albumin present in the well during the immunological incubation.

The principle of the assay is summarised below:



MATERIALS AND EQUIPMENT REQUIRED

In addition to standard laboratory equipment, the following material is required:

- ☞ Precision micropipettes (20 to 1000 μL)
- ☞ Spectrophotometer plate reader (405 or 414 nm filter)
- ☞ Microplate washer (or washbottles)
- ☞ Microplate shaker
- ☞ Distilled or deionized water
- ☞ Polypropylene tubes

SAMPLE COLLECTION & PREPARATION

This assay may be used to measure albumin in rat urine sample. Urine samples are collected in propylene tubes. Urine are collected and kept at -20°C until assay.

No prior extraction procedure is necessary to measure albumin in urine samples.



Explore our innovative technologies for **YOUR** research

REAGENT PREPARATION

The coated plates and reagents are provided ready to use.

☞ EIA buffer

Reconstitute one vial with 50 mL of distilled or deionized water. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion. Stability at 4°C: 1 month.

☞ Rat albumin standard

Reconstitute the vial with 1 mL of distilled or deionized water. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion. The concentration of the first standard is 1000 ng/mL. Prepare seven propylene tubes (for the seven other standards) and add 500 μ L of EIA buffer into each tube. Add 500 μ L of the first tube (containing the first standard) to the second tube. Continue this procedure for the other tubes. Thus, standard concentrations are: 1000 (S1), 500 (S2), 250 (S3), 125 (S4), 62.5 (S5), 31.25 (S6), 15.63 (S7) and 7.81 ng/mL (S8), respectively. Stability at 4°C: 1 week.

☞ Quality Control

Reconstitute one vial with 1 mL of distilled or deionized water. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion. Stability at 4°C: 1 week.

☞ Rat albumin-AChE tracer

Reconstitute one vial with 5 mL of EIA buffer. Allow it to stand 5 minutes until completely dissolved and then mix thoroughly by gentle inversion. Stability at 4°C: 1 month.

☞ Rat albumin antiserum

Reconstitute one vial with 5 mL of EIA buffer. Allow it to stand 5 minutes completely dissolved and then mix thoroughly by gentle inversion. Stability at 4°C: 1 week.

☞ Wash buffer

Dilute 1 mL of concentrated Wash buffer to 400 mL with distilled or deionized water. Add 200 μ L of tween 20 (Use a magnetic stirrer to mix the contents). Stability at 4°C: 1 week.

☞ Ellman's Reagent

Five minutes before use, reconstitute with 50 mL of distilled or deionized water. The tube contents should be thoroughly mixed. Stability at 4°C and in the dark: 4 days.

ASSAY PROCEDURE

It is recommended to perform the assays in duplicate and to follow the instructions hereafter.

PLATE PREPARATION

Prepare the wash buffer as indicated in the reagent preparation section. Open the plate packet and select the sufficient strips for your assay and place the unused strips back in the packet (stored at 4°C). Rinse each well five times with the wash buffer (300 μ L/well).

Just before distributing reagents and samples, remove the buffer from the wells by inverting the plate and shaking out the last drops.

DISTRIBUTION OF REAGENTS AND SAMPLES

A plate set-up is suggested on the following page. The contents of each well may be recorded on the sheet provided with the kit.

PIPETTING THE REAGENTS

Note that the first column should be left empty for blanking Ellman's reagent. All samples and reagents must reach room temperature prior to performing the assay. Use different tips to pipet the buffer, standard, sample, tracer, antiserum and other reagents.

	12	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	11	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	10	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	9	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	8	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	7	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	6	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	5	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)	(*)
	4	(S5)	(S5)	(S6)	(S6)	(S7)	(S7)	(S8)	(S8)	(S8)	(S8)	(S8)
	3	(S1)	(S1)	(S2)	(S2)	(S3)	(S3)	(S4)	(S4)	(S4)	(S4)	(S4)
	2	(NSB)	(NSB)	(Bo)	(Bo)	(Bo)	(Bo)	(*)	(*)	(*)	(*)	(*)
	1	(B)	(B)	(B)	(B)	(B)	(B)	(B)	(B)	(B)	(B)	(B)
A			B	C	D	E	F	G	H			

B : Blank
 NSB : Non Specific Binding
 Bo : Maximum Binding
 S1-S8 : Standards 1-8
 * : Samples or Quality Controls

- ↪ EIA buffer: Dispense 100 µL to Non Specific Binding (NSB) wells and 50 µL to Maximum Binding (Bo) wells.
- ↪ Rat albumin standard: Dispense 50 µL of each of the eight standards (S1 to S8) in duplicate to appropriate wells. Start with the lowest concentration standard (S8) and equilibrate the tip in the next higher standard before pipetting.
- ↪ Quality Control and samples: Dispense 50 µL in duplicate to appropriate wells. Highly concentrated samples may be diluted in EIA buffer.
- ↪ Rat albumin AChE tracer: Dispense 50 µL to each well except the Blank (B) wells.
- ↪ Rat albumin antiserum: Dispense 50 µL to each well except the Blank (B) wells and the Non Specific Binding (NSB) wells.

INCUBATING THE PLATES

Cover the plate with a plastic film and incubate for 16-20 hours at 4 °C.

DEVELOPING AND READING THE PLATE

Reconstitute the wash buffer and Ellman's Reagent as indicated in reagent preparation section. Empty the plate by turning over and shaking. Then, wash each well five times with the wash buffer (300 µL/well). Dispense 200 µL of Ellman's Reagent to the 96 wells. Incubate in the dark (plate covered with an aluminium sheet) at room temperature. Optimal development is obtained using an orbital shaker. The plate should be read between 405 and 414 nm (yellow colour) when the Maximum Binding (Bo) wells reach an absorbance of 0.2-0.8 unit.



Explore our innovative technologies for YOUR research

Enzyme Immunoassay Protocol (Volumes are in μL)					
	Blank	Non specific binding	Maximum binding	Standard	Sample
Buffer	-	100	50	-	-
Standard	-	-	-	50	-
Sample	-	-	-	-	50
Tracer	-	50	50	50	50
Antiserum	-	-	50	50	50
Cover the plate, incubate at 4 °C for 16-20h					
Wash the plate 5 times					
Ellman's reagent	200	200	200	200	200
Incubate the plate with an orbital shaker in the dark at room temperature					
Read the plate between 405 and 414 nm					

DATA ANALYSIS

Make sure that your plate reader has subtracted the absorbance readings of the Blank well (absorbance of Ellman's reagent) from the absorbance readings of the rest of the plate. If not, do it now.

- ↪ Calculate the average absorbance for each NSB, Bo, standard and sample.
- ↪ Calculate the B/Bo (%) for each standard and sample: (average absorbance of standards or sample - average absorbance of NSB) divided by (average absorbance of Bo - average absorbance of NSB) & multiplied by 100.
- ↪ Using a semi-log graph paper, plot the B/Bo (%) for each standard point (y axis) versus the concentration (x axis). Draw a best-fit line through the points.
- ↪ To determine the concentration of your samples, find the B/Bo (%) value on the y axis. Read the corresponding value on the x axis which is the concentration of your unknown sample. Samples with a concentration superior to 500 ng/mL should be re-assayed after dilution in EIA buffer.
- ↪ Most plate readers are supplied with curve-fitting software capable of graphing this type of data (logit/log or 4-parameter). If you have this type of software, we recommend using it. Refer to it for further information.



Explore our innovative technologies for **YOUR** research

TYPICAL DATA

EXAMPLE DATA


The following data are for demonstration purpose only. Your data may be different and still correct. These data were obtained using all reagents as supplied in this kit under the following conditions: 1 hour developing at 20°C, reading at 414 nm. A spline fitting was used to determine the concentrations.

	mAU	B/Bo (%)
NSB	1	
Bo	520	100
Standard 1000 ng/mL	34	6.36
Standard 500 ng/mL	59	11.2
Standard 250 ng/mL	101	19.3
Standard 125 ng/mL	154	29.5
Standard 62.5 ng/mL	226	43.4
Standard 31.25 ng/mL	292	56.1
Standard 15.63 ng/mL	385	74.0
Standard 7.81 ng/mL	437	84.0
QC	279	

ACCEPTABLE RANGE

- ☞ Bo absorbance: > 200 mAU in the conditions indicated above.
- ☞ Ratio NSB absorbance / Bo absorbance: < 0.1.
- ☞ 50% B/Bo (%): 22 to 42 ng/mL (mean: 32 ng/mL).
- ☞ QC sample: see the label on the vial.

RAT ALBUMIN STANDARD CURVE





Explore our innovative technologies for YOUR research

ASSAY VALIDATION AND CHARACTERISTICS

The Enzyme Immunoassay of rat albumin has been fully validated for its use in rat urine. Assay characteristics are the following :

- ☞ Cross-reactivity with serum albumin bovine: <0.01%.
- ☞ The limit of detection calculated as the concentration of albumin corresponding to the NSB average minus three standard deviations: 13 ng/mL.
- ☞ Intra-assay variation on 4 QCs (n=8):
 - 300 ng/mL: 1.5 %
 - 150 ng/mL: 7.1 %
 - 75 ng/mL: 4.6 %
 - 30 ng/mL: 4.2 %
- ☞ Inter-assay variation on 4 QCs (n=8):
 - 300 ng/mL: 4.8 %
 - 150 ng/mL: 4.1 %
 - 75 ng/mL: 3.8 %
 - 30 ng/mL: 8.2 %
- ☞ Accuracy on 4 QCs (n=8):
 - 300 ng/mL: 94 %
 - 150 ng/mL: 91 %
 - 75 ng/mL: 98 %
 - 30 ng/mL: 91 %
- ☞ Limit of quantification: 20 ng/mL

ASSAY TROUBLE SHOOTING

- ☞ Bo value is too low: incubation in wrong conditions (time or temperature) or reading time too short or Rat albumin-AChE tracer, Rat albumin antiserum or Ellman's reagent have not been dispensed.
- ☞ NSB value too high: contamination of NSB wells with Rat albumin antiserum or inefficient washing.
- ☞ High dispersion of duplicates: poor pipetting technique or irregular plate washing.
- ☞ IC₅₀ or QC concentrations not within the expected range: wrong preparation of standards.
- ☞ Analyses of two dilutions of a biological sample do not agree: interfering substances are present. Sample must be purified prior to EIA analysis.

These are a few examples of trouble shooting that may occur. If you need further explanation, SPI-BIO will be happy to answer any questions or information about this assay. Please feel free to contact our technical support staff by letter, phone (33 (0)1 39 30 62 60), fax (33 (0)1 39 30 62 99) or E-mail (sales@spibio.com), and be sure to indicate the lot number of the kit (see outside of the box).

SPI-BIO offers a training workshop in EIA practice & theory. This workshop is given twice a year. For further information, please contact our Customer Relation Representative (33 (0)1 39 30 62 60).

BIBLIOGRAPHY

- ☞ **Grassi J. & Pradelles Ph.**
Compounds labelled by the acetylcholinesterase of *Electrophorus Electricus*. Its preparation process and its use as a tracer or marker in enzyme-immunological determinations.
United States patent, N°1,047,330. September 10, 1991
- ☞ **Vienet R., Grognet J.-M., Ezan E., Lecaque D., Hamon G. & Corman B.**
Effect of chronic converting-enzyme inhibition on kidney function of senescent hypertensive rats.
Journal of cardiovascular pharmacology Vol. 23 (suppl. 4), 1994, S19-S25